

THE GOVERNMENT PHARMACEUTICAL ORGANIZATION

RAW MATERIAL SPECIFICATION

Title: Hydrogen Peroxide Solution [50% w/w]

Spec. No. : SP-AK30-H3

(Item No. 41032100)

Reference(s): GPO Specification

Rev. No. : 04

Other Requirements: BP 2020 page I-1265 to I-1266

Page

: 1/2

[Hydrogen Peroxide Solution 30% w/w]

GPO Specification

Test Parameters	Requirement			
Description	Colourless, clear liquid.			
Identification	A. The solution becomes colourless with evolution of gas.			
	B. To 1 mL, add 0.1 mL of dilute hydrochloric acid R and 0.1 mL of potassium iodide solution.			
	A brown colour appears. Black particles may be formed.			
Acidity	To 10 mL, add 100 mL of purified water and 0.25 mL of methyl red solution. Not less than			
	0.05 mL and not more than 0.5 mL of 0.1 M sodium hydroxide is required to change the			
	colour of the indicator.			
Organic stabilisers	The residue weighs not more than 10 mg.			
Non-volatile residue	The residue weighs not more than 20 mg.			
Assay	48.0% w/w to 51.0% w/w of Hydrogen peroxide (H ₂ O ₂).			

เอกสารไม่ควบกุม ใช้ในการจัดซื้อ

Prepared by:	Reviewed by:		Approved by:	Eff. Date
	Tournes Lumbront, 14100120	// 18/06/20	Victoria Ruemanaroi, 18106/20	15/07/20
J. Noppadech 17/01/20			Caangreen 103 / 1915 125	13/07/20
Head of Raw Material	Director of Raw Material	Director of Drug Registration and	Director of Quality Assurance	
Standard Section 2	Standard Division	Pharmacovigilance Division	Department (Alive)	

LY respondent CUTY NO.

THE GOVERNMENT PHARMACEUTICAL ORGANIZATION

RAW MATERIAL SPECIFICATION

Title: Hydrogen Peroxide Solution [50% w/w]

Spec. No. : SP-AK30-H3

(Item No. 41032100)

Reference(s): GPO Specification

Rev. No. : 04

Other Requirements: BP 2020 page I-1265 to I-1266

Page

: 2/2

[Hydrogen Peroxide Solution 30% w/w]

Product Information

Approved source(s)	Refer to current version of Approved Supplier List of Hydrogen Peroxide Solution [50% w/w] (Item No. 41032100).		
Sampling plan	For Identification: 100%. For Other Tests: n plan.		
Testing procedure	Tests to be performed as per current version of WI-AK30-H3.		
Storage condition	Store at a condition stated on the label from the manufacturer. Protected from light, and if the solution does <u>not</u> contain a stabiliser, at a temperature below 15 °C (It should not be stored for long period).		
Retest period	1 year after first testing date.		

History of changes

Rev. No.	Description	Effective Date
02	ทบทวนเอกสาร อ้างอิงตาม BP2011 และ GPO Specification	19/08/11
03	ทบทวนเอกสาร โดยอ้างอิงตาม GPO Specification และ BP 2014 [Hydrogen Peroxide	25/10/16
	Solution 30% w/w]	
04	ทบทวนเอกสาร โดยอ้างอิงตาม GPO Specification และ BP 2020 [Hydrogen Peroxide	15/07/20
	Solution 30% w/w]	
	เอกสารในควบคุม	

ใช้ในการจัดซื้อ

Prepared by:	Reviewed by:		Approved by:	Eff. Date
	tarnia	1-11	lichim	
J. Noppodech/17/01/20	Lumbount, 14101120	18/06/20	Rungunga, 18106/20	15/07/20
Head of Raw Material	Director of Raw Material	Director of Drug Registration and	Director of Quality Assurance	2 25 45
Standard Section 2	Standard Division	Pharmacovigilance Division	Department (Auning)	